


GAS REGULATORS




Click-Style Regulator with Yoke



Click-Style Regulator with Nut & Nipple



Compact Regulator with Yoke



Compact Regulator with Nut & Nipple

Amvex Gas Regulators are available in two different styles: 50 PSI Preset Compact Regulator and 50 PSI Preset Click Style Regulator. They are offered in a variety of gases and come with either a CGA Nut and Nipple or CGA Yoke cylinder connection. The Click Style Regulator has a built in flowmeter and the Compact Regulator can have a tube style flowmeter added if required.

The Regulators are designed to meet all current CGA, HCA, FDA and CE requirements. They are cost effective, durable and easy to use. They are all 100% tested and quality inspected prior to shipping. All regulators come with a 5 year warranty.

Available in a variety of colors



FEATURES AND BENEFITS

Click Style Regulator

- Amvex Click Style Regulators are available in 50 PSI Preset
- Available with a CGA Nut and Nipple or CGA Yoke Style
- Available in Oxygen only
- Aluminum body with brass core to ensure safety
- Built in Flowmeter available in 0-4 LPM, 0-8 LPM, 0-15 LPM, and 0-25 LPM
- Available in a variety of colors (blue is standard)



Compact Regulator

- Amvex Compact regulators are available in 50 PSI Preset
- Available with a CGA Nut and Nipple or CGA Yoke Style
- Available in Oxygen, Air, N2O, and CO2
- Available with flowmeter options
- Machined brass
- Available in USA and ISO colors

PART NUMBERS

Click Style Regulators

Nut & Nipple

CGA 540 25 LPM CLICK: **GR25-540**
CGA 540 15 LPM CLICK: **GR15-540**
CGA 540 8 LPM CLICK: **GR8-540**
CGA 540 4 LPM CLICK: **GR4-540**

Yoke Style

CGA 870 25 LPM CLICK: **GR25-870**
CGA 870 15 LPM CLICK: **GR15-870**
CGA 870 8 LPM CLICK: **GR8-870**
CGA 870 4 LPM CLICK: **GR4-870**

LPM Increments

25 LPM: 0.5,1,2,3,5,6,8,10,15,25
15 LPM: .12,.5,1,1.5,2,3,4,6,8,10,15
08 LPM: .12,.5,1,1.5,2,2.5,3,4,5,6,8
04 LPM: .03,.06,.12,.25,.5,.75,1,1.5,2,3,4

Click Style regulator color options available are: Purple "P", Green "G", Gold "Y", Red "R", Black "B". Please add color letter after the part number to specify color option. Example for color part is: **GR15-870P**.

Compact Regulators

Nut & Nipple and Yoke Style

CGA 540 (OXY NUT & NIPPLE): **GR-540**
CGA 870 (OXY YOKE STYLE): **GR-870**
CGA 346 (AIR NUT & NIPPLE): **GR-346**
CGA 950 (AIR YOKE STYLE): **GR-950**
CGA 326 (N2O NUT & NIPPLE): **GR-326**
CGA 910 (N2O YOKE STYLE): **GR-910**
CGA 320 (CO2 NUT & NIPPLE): **GR-320**
CGA 940 (CO2 YOKE STYLE): **GR-940**

With Flowmeters

CGA 540 (WITH 15 LPM FLOWMETER): **GR-540FM**
CGA 870 (WITH 15 LPM FLOWMETER): **GR-870FM**
CGA 346 (WITH 15 LPM FLOWMETER): **GR-346FM**
CGA 950 (WITH 15 LPM FLOWMETER): **GR-950FM**
CGA 326 (WITH 15 LPM FLOWMETER): **GR-326FM**
CGA 910 (WITH 15 LPM FLOWMETER): **GR-910FM**
CGA 320 (WITH 15 LPM FLOWMETER): **GR-320FM**
CGA 940 (WITH 15 LPM FLOWMETER): **GR-940FM**

