


AMICO ARTWALL GAS SYSTEM (M/N:W-ARTWALL-03)



TYPE:

QUANTITY:

SYMBOL	QTY.	DESCRIPTION
SL	1	VACUUM SLIDE WALL MOUNTED
O	1	SWIVEL TYPE GAS OUTLET, OXYGEN
V	1	SWIVEL TYPE GAS OUTLET, VACUUM
A	1	SWIVEL TYPE GAS OUTLET, AIR
ASB	1	SPHYMOMANOMETER C/W HARDWARE
		ANERIOD (OPTIONAL)

FEATURES:

- RECESSED ENCLOSURE (Ga#16 Galvanized steel)
- SWIVEL TYPE GAS OUTLETS
- FLOWMETERS, REGULATORS CAN BE ATTACHED TO THE SWIVEL OUTLETS.
- HORIZONTAL MOVEMENT OF PICTURE FRAME .
- ACCEPTS FRAMED PRINTS.

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING INFORMATION.

APPROVAL SIGNATURE _____

DATE _____

PHONE NO. _____

AMICO Corporation
 85 Fulton Way
 Richmond Hill, Ontario
 L4B 2N4, CANADA
 Toll-Free: 1-877-462-6426
 Tel: (905) 764-0800
 Fax: (905) 764-0862
 www.amico.com

HOSPITAL _____

TYPICAL DRAWING

LOCATION _____

QTY. (_____) TYPE A UNITS AS SHOWN / (_____) TYPE B UNITS OPPOSITE

A. NURSE CALL MFGR: _____ MODEL NO.: _____

B. MEDICAL GAS MFGR.: AMICO TYPE CONNECTION: _____

C. FINISH: _____ CEILING HEIGHT: _____

DRWG. NO. _____

YYMMDD-### _____

DRAWN BY: _____

CHECKED BY: CS _____

REV. NO.: _____

DATE: MM/DD/YY _____