


# AMICO ARTWALL GAS SYSTEM (M/N:W-ARTWALL-05)



TYPE:

QUANTITY:

SYMBOL	QTY.	DESCRIPTION
SL	2	VACUUM SLIDE WALL MOUNTED
O	2	SWIVEL TYPE GAS OUTLET, OXYGEN
V	2	SWIVEL TYPE GAS OUTLET, VACUUM
A	1	SWIVEL TYPE GAS OUTLET, AIR
ASB	1	SPHYMOMANOMETER C/W HARDWARE
		ANEROID (OPTIONAL)

FEATURES:

- RECESSED ENCLOSURE (Ga#16 Galvanized steel)
- SWIVEL TYPE GAS OUTLETS
- FLOWMETERS, REGULATORS CAN BE ATTACHED TO THE SWIVEL OUTLETS.
- HORIZONTAL MOVEMENT OF PICTURE FRAME .
- ACCEPTS FRAMED PRINTS.

**IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING INFORMATION.**

APPROVAL SIGNATURE \_\_\_\_\_

DATE \_\_\_\_\_

PHONE NO. \_\_\_\_\_



85 Fulton Way  
Richmond Hill, Ontario  
L4B 2N4, CANADA  
Toll-Free: 1-877-462-6426  
Tel: (905) 764-0800  
Fax: (905) 764-0862

www.amico.com

HOSPITAL		<b>TYPICAL DRAWING</b>	
LOCATION		A. NURSE CALL MFGR: _____ MODEL NO.: _____	
QTY.	(____) TYPE <u>A</u> UNITS AS SHOWN / (____) TYPE <u>B</u> UNITS OPPOSITE	B. MEDICAL GAS MFGR.: <u>AMICO</u> TYPE CONNECTION: _____	
		C. FINISH: _____ CEILING HEIGHT: _____	

DRWG. NO.
YYMMDD-###
DRAWN BY:
CHECKED BY:CS
REV.NO.:
DATE:MM/DD/YY