



AMICO FLATWALL - 19" Unit (M/N: FW1907R-GEN-P)
 THREE QUARTER LENGTH GEN. CARE SEMI RECESSED



TYPE:
 QUANTITY:
 FLATWALL 19"

SYMBOL	QTY.	DESCRIPTION
FW	1	BASIC FLATWALL 19" WIDTH
N2	2	RECEPTACLE - DUPLEX IVORY
E2	3	RECEPTACLE - DUPLEX RED
NC	1	NURSE CALL PROVISION
CB	1	CODE BLUE PROVISION
SLM	1	VACUUM SLIDE WITH MOUNT
A	1	GAS, AIR
O	1	GAS, OXYGEN
V	1	GAS, VACUUM

IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATIONS ARE CORRECT, AND SUPPLY THE REMAINING INFORMATIONS.

APPROVAL SIGNATURE _____

DATE _____

PHONE NO. _____

HOSPITAL	85 Fulton Way Richmond Hill, Ontario L4B 2N4, CANADA Tel: (905) 1-877-462-6426 Fax: (905) 764-0800 www.amico.com	DRWG. NO. 031120-C02
LOCATION		DRAWN BY: CE
QTY: (A) TYPE _____ UNITS AS SHOWN / (B) TYPE _____ UNITS OPPOSITE		CHECKED BY: CS
A. NURSE CALL MFGR.: _____	MODEL NO.: _____	DATE: 11/20/03
B. MEDICAL GAS MFGR.: AMICO	TYPE CONNECTION: _____	
C. FINISH: _____	CEILING HEIGHT: _____	

