

HORIZONTAL HEADWALL



MAJESTIC SERIES

A Corporation
AMICO


Majestic Series Horizontal Headwall System

Amico's Majestic Series of Horizontal Headwalls is an attractive, easy to maintain solution for providing medical services to various acuity levels. The Majestic Headwall system simplifies equipment management at bedside while de-institutionalizing the look of essential services within the care area. Available in single and multiple tier layouts, the Majestic series headwall system is built to suit the end user's needs.



- An integrated system, the Majestic series headwall system also includes overbed lights, Bed Locators and a full line of equipment management accessories.
- Custom designed to suit the needs of your facility - almost any size and configuration of outlets can be accommodated.
- Factory wired and piped to facilitate single point site connection, thus compressing traditional construction schedules.
- Great for new construction or renovations. Designed for back wall or chase feeding to suit the needs of your project.
- Built in electrical race ways for ease of wiring and future expansion.

Amico Overbed Light




Amico Bed Locator

Standard color is Snowflake White. Other colors available upon request.

Easy to clean, vacuum formed Kydex-T surface.

Accessories



I.V. Hook
Used to hold I.V. solution bags. Duplex and Quad Hooks Available.

W-A-IVHOOK2-H



Aneroid Manometer
Aneroid with large easy to read dial makes blood pressure reading convenient and accurate. Cuff basket available. Mercury free.

W-A-AMANOM-CB-H



Channel Bracket
It's used for mounting various 1/2" diameter post accessories, such as: I.V. Hooks, Exam Lights and Swivel Shelves.

W-A-HTA-05




Vacuum Bottle Slide
Used to hold suction canisters and other secondary equipment.

W-A-SLIDE-H



Chart Holder
Keeps patient information within reach.

W-A-CHARTHLD-H



Storage Basket
Convenient storage of medical supplies and staff protection products. Wide variety of sizes are available.

W-A-BB-090910-H (for this size and mount)



Universal Sharps Holder
Holds various shapes and sizes of sharp containers.

W-A-SPHOLDER-H



Resuscitation Bag Holder
Holds a resuscitation bag within easy reach for the medical staff.

W-A-RESUHOLD-H



Tubing Clips
Provides access and organization of the bedside tubing and cables. 1 to 4 clips are available.

W-A-CLIPTUB-H-2 (for this size and mount)




Swivel Shelf
Used for writing or supply storage. It can hold up to 25lbs.


W-A-SHELS-1212H

* For a complete listing of accessories, please refer to our Accessory Brochure.


Headwall Finishes



Antique White Powder Coated



Clear Anodized Aluminum




Champagne Anodized Aluminum



Light Bronze Anodized Aluminum


Other powder coated finishes and anodized aluminum finishes are available upon request. Due to printing limitations, colors may vary slightly. Please contact Amico for color samples.

DISTRIBUTED BY:



Laminate Inserts

Wide variety of laminate available.



Amico products are listed to both Canadian and U.S. requirements.



85 Fulton Way, Richmond Hill, ON L4B 2N4 Canada • Tel: 1-877-462-6426 • Fax: 905-764-0862
14 Madison Road, Fairfield, NJ 07004 USA • sales@amico.com • www.amico.com