

**TABLE FOR CONVERSION PRODUCTS**

Example: O-FCO1-QD-L-GAS → O-FCO1-QD-U-OXY

For the "L", please specify the language:

- U = English (NFPA) Oxy, Air, Vac and WAG only
- E = English (CSA)
- F = French (CSA)
- S = Spanish (NFPA)

For the "GAS" type required, please specify from the following:

- Oxygen = OXY
- Medical Air = AIR
- Nitrous Oxide = N2O
- Carbon Dioxide = CO2 \* Only Available in DISS or QD
- Nitrogen = NIT \*\* Only Available in DISS
- MedVac = VAC
- Waste Anesthetic Gas Disposal = WAG (NFPA)
- Anesthetic Gas Scavenging System = AGS (CSA)

		OHIO MEDICAL/OHMEDA				PURITAN BENNETT			NCG	OXEQUIP/MEDSTAR	
		DIAMOND 1	DIAMOND 2		DIAMOND 3	SERIES 73	SERIES "A"	LITHARGE	SERIES 148,400,460	SERIES # 7 OR 5/8-24	SMOOTH BORE OR 5/8-27
			Removing rough-in parts	Using existing rough-in parts							
<b>FRONT ASSEMBLY</b>											
For each outlet choose one of:	<b>QD</b>	O-FCO1-QD-L-GAS	O-FCO1-QD-L-GAS	O-FCO3-QD-L-GAS	O-FCO3-QD-L-GAS	O-CPBF-QD-L-GAS	O-CPBF-QD-L-GAS	O-CPBF-QD-L-GAS	O-FASC-QD-L-GAS	O-FASC-QD-L-GAS	O-CPBF-QD-L-GAS
	<b>DI</b>	O-FCO1-DI-L-GAS	O-FCO1-DI-L-GAS	O-FCO3-DI-L-GAS	O-FCO3-DI-L-GAS	O-CPBF-DI-L-GAS	O-CPBF-DI-L-GAS	O-CPBF-DI-L-GAS	O-FASC-DI-L-GAS	O-FASC-DI-L-GAS	O-CPBF-DI-L-GAS
	<b>CH</b>	Not Available	Not Available	Not Available	Not Available	Not Available	Not Available	Not Available	O-FASC-CH-L-GAS	O-FASC-CH-L-GAS	Not Available
	<b>PB</b>	Not Available	Not Available	Not Available	Not Available	O-CPBF-PB-L-GAS	O-CPBF-PB-L-GAS	O-CPBF-PB-L-GAS	O-FASC-PB-L-GAS	O-FASC-PB-L-GAS	O-CPBF-PB-L-GAS
<b>BACK ASSEMBLY</b>		O-VC-OH1-L-GAS	O-VC-OH2-GAS	Not needed	Not needed	O-CPBBAK-L-GAS	O-CPBBAK-L-GAS	O-CPBBAK-L-GAS	O-VC-NCG-L-GAS	O-VC-OX1-L-GAS	O-MS-SBCO-L-GAS
<b>CONNECTOR</b>		Not needed	Not needed	Not needed	Not needed	O-CPBCON-73	O-CPBCON-SA	O-CPBCON-LI	Not needed	Not needed	Not needed
<b>SPACERS</b>						O-CON-FIL-PB325 or O-CON-FIL-PB500			O-CON-FIL-PB500	INCLUDED IN BACK ASSEMBLY	
Single		O-CON-FIL-OH1-1	Not needed	Not needed	Not needed						
Double		O-CON-FIL-OH1-2	Not needed	Not needed	Not needed						
Triple		O-CON-FIL-OH1-3	Not needed	Not needed	Not needed						
		Above spacers are painted white.				Above spacers are chrome finish.					
<b>FINISH FRAME</b>		O-TRIM-M325C	O-TRIM-M450C	O-TRIM-M450C	O-TRIM-M450C	O-TRIM-M325C or O-TRIM-M500C			O-TRIM-M500C	INCLUDED IN BACK ASSEMBLY	
H" x W"		(5.5"x3.25")	(4.5"x4.5")	(4.5"x4.5")	(4.5"x4.5")	(5.5"x3.25") or (5.5"x 5.0")			(5.5"x 5.0")		
<b>INSTALLATION TOOL</b>		O-VC-OH1-TOOL	Not needed	Not needed	Not needed	O-CPBTOOLKIT	Not needed	Not needed	O-VC-OH1-TOOL	O-VC-OH1-TOOL	Not needed