



AMICO RAIL ASSEMBLY

SPECIFICATIONS:

Wall Mounted Amico Rail System

- 2.0' (609) W-RAIL-A-024
- 3.0' (914) W-RAIL-A-036
- 4.0' (1219) W-RAIL-A-048
- 5.0' (1524) W-RAIL-A-060
- 6.0' (1829) W-RAIL-A-072
- 8.0' (2438) W-RAIL-A-096
- 10.0' (3048) W-RAIL-A-120

The Amico Rail System is a complete wall mounted equipment storage system. Installed at each patient location, it provides ready access to commonly used equipment such as suction regulators, flowmeters, exam lights, shelves and collection bottles.

The Amico Rail System accommodates a wide variety of snap-in brackets which interface with various equipment. The ARS is permanently attached to the wall in the vicinity of the Medical Gas Outlets. The smooth unobstructed design of the ARS's lower surface facilitates easy cleaning.

Each ARS includes an aluminum mounting bar which is field attached to the building structure by the customer. The ARS is available in a variety of lengths and can be field modified to job-site conditions if required.



PROJECT: _____

D I S T R I B U T E D B Y :

A Corporation
AMICO 14 Madison Rd., Fairfield, NJ 07004 USA
85 Fulton Way, Richmond Hill, ON L4B 2N4 Canada

Tel: 1-877-GOAMICO (1-877-462-6426)
Fax: (905) 764-0862 www.amico.com