


# AMICO HORIZONTAL MEDICAL WALL SYSTEM MAJESTIC SERIES 1 SINGLE TIER WITH NO CHASE



TYPE:  
QUANTITY:

SYMBOL	QTY.	DESCRIPTION
HC	2	HORIZONTAL CHANNEL
N2	2	DUPLEX RECEPTACLE - IVORY
E2	2	DUPLEX RECEPTACLE - RED
NC	1	NURSE CALL PROVISION
CB	1	CODE BLUE PROVISION
DA	1	DATA PROVISION
O	2	GAS, OXYGEN
A	2	GAS, AIR
V	3	GAS, VACUUM
SLM	3	VACUUM SLIDE C/W MOUNT

\*BACK BOX SIZE (NC):

BED LOCATOR SYSTEM DETAILS		
DEEP	WIDE	HIGH
SYMBOL	QTY.	DESCRIPTION
N2	2	DUPLEX RECEPTACLE - IVORY
E2	2	DUPLEX RECEPTACLE - RED
TL	1	TELEPHONE PROVISION

**IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION ARE CORRECT, AND SUPPLY THE REMAINING INFORMATION.**

APPROVAL SIGNATURE \_\_\_\_\_ DATE \_\_\_\_\_ PHONE NO. \_\_\_\_\_



HOSPITAL	TYPICAL DRAWING
LOCATION	
QTY.	(_____) TYPE _____ UNITS AS SHOWN / (_____) TYPE _____ UNITS OPPOSITE

A. NURSE CALL MFGR: _____ MODEL NO.: _____
B. MEDICAL GAS MFGR.: _____ AMICO TYPE CONNECTION: _____
C. FINISH: _____ CEILING HEIGHT: _____

DRWG. NO.	031104-C03
DRAWN BY:	CE
CHECKED BY:	CS
REV. NO.:	00
DATE:	11/04/03